

Name _____

Date _____

1 What test results are included and reported for the Alere Afinion ACR assay?

- a. ACR
- b. ACR and Albumin
- c. Creatine
- d. ACR, Albumin, and Creatinine

2 What are acceptable sample types for the Alere Afinion ACR assay?

- a. Fresh Urine
- b. Frozen Urine
- c. Serum
- d. a. and b.

3 How long are Alere Afinion ACR controls stable for once the vials are opened?

- a. 6 weeks
- b. 8 weeks
- c. 12 weeks
- d. Until the expiration date printed on the box and vials

4 What does an ACR value of “---” mean?

- a. Patient has kidney failure
- b. ACR value is too high to be reported
- c. Albumin and/or Creatinine value is outside of measureable range
- d. Alere Afinion Analyzer needs service call technical support

5 What is the sample size for the Alere Afinion ACR assay?

- a. 1.5 uL
- b. 2.5 uL
- c. 3.5 uL
- d. 5.0 uL

6 True or False: Alere Afinion ACR Controls can be stored at room temperature (59-77°F)

7 Urine samples can be stored refrigerated for how many days?

- a. 1 day
- b. 5 days
- c. 10 days
- d. Samples cannot be stored refrigerated

8 How often should Alere Afinion ACR Controls be run as per the manufacturer's recommendations?

- a. At least every 30 days
- b. With each new shipment of ACR test kits
- c. With each new lot of ACR test kits
- d. All of the above

9 True or False: Alere Afinion ACR test cartridges can be used immediately after removing them from the refrigerator?

10 How long are Alere Afinion ACR test cartridges stable for if stored at room temperature (59-86°F)?

- a. 5 hours
- b. 8 hours
- c. 24 hours
- d. Until the expiration date printed on the box and foil pouch



Knowing now matters.™

Alere Technologies AS
Kjelsåsveien 161
P.O. Box 6863 Rodeløkka
NO-0504 Oslo, Norway