

BinaxNOW® Influenza A & B Quiz

Name: _____

Date: _____

Circle T (True) or F (False) for each Question:

1. The BinaxNOW Influenza A& B test kit requires refrigeration at 2° to 8°C. T F
2. The BinaxNOW Influenza A & B test device should not be removed from the foil pouch until just before use. T F
3. Throat swabs are an acceptable sample type for testing on the BinaxNOW Influenza A & B test. T F
4. Nasal Wash samples may be stored at room temperature for 24 hours. T F
5. Transport media may be used with the BinaxNOW Influenza A & B test. T F
6. NP swabs may be stored at room temperature for 24 hours before elution. T F
7. Sample should be added quickly to the top of the white sample pad. T F
8. Test Results should be read at 15 minutes. T F
9. If the Control Line remains BLUE in color the test results are valid. T F
10. The appearance of a pink-to-purple Control Line and a pink-to-purple Sample Line is a positive result. T F

Reviewed By: _____

Date: _____